

# Nu-Tek Yeast Peptone SY6B

Product Code: 2200

## Maximum Microbial Growth with Short Production Cycles

Peptone SY6B is a carefully balanced non-animal peptone developed for maximum microbial growth with short production cycles. Peptone SY6B was designed to be efficient for fermentation and diagnostic testing. It is especially effective when used for growth of *E. coli*, *Bifidobacterium sp.* and *Lactobacillus sp.* When used at typical use levels of yeast extracts and peptones in microbial fermentations, SY6B will generally reach optimum cell growth yields in a shorter cycle time with comparable or greater yields.

Due to its cost effectiveness and excellent uniformity, Peptone SY6B can be used at all levels of product development including research, pilot plant, seed fermentation and commercial fermentation. In addition it is animal-free and kosher certified.

### SPECIFICATIONS:

#### Typical Analysis:

Moisture	< 5.0%
Total Nitrogen (dry basis)	8.0%–9.5%
Amino Nitrogen (dry basis)	3.2%–4.8%
AN/TN Ratio	0.34–0.6
Total Carbohydrates	30.0%–34.0%
Ash	< 20.0%
pH	5.0–6.5
Clarity	< 20 NTU

#### Microbiological Analysis:

SPC	< 5,000 cfu/g
Coliform	< 5 cfu/g
<i>E. coli</i>	< 0.3 mpn/g
Yeast & Mold	< 50 cfu/g
Salmonella	Negative/25g

#### Packaging:

25 Kg polyethylene bag in a corrugated box.

#### Storage Requirements:

This product should be kept in a cool, dry and ventilated place. The product should be used as soon as possible after the bag has been opened. Any partially used bags must be sealed properly to prevent moisture adsorption by product.

#### Shelf Life:

The shelf life of the product is 3 years when stored at < 30°C.



*This information is presented in good faith and great care was used in its preparation. However, no warranty, guarantee, or freedom from patent infringement is implied or inferred. This information is offered solely for your consideration and verification.*

technologies for the biotech and nutraceutical ingredient industries

952.936.3600

5400 Opportunity Court,  
Suite 120  
Minnetonka, MN 55343

[www.Nu-TekBioScience.com](http://www.Nu-TekBioScience.com)